

SOP-DRYING-001: Fluid-Bed Drying Step (Synthetic)

1. Purpose. This SOP defines the drying step parameters for oral solid dosage batches processed in Oven-01 through Oven-06.
2. Process parameters. Target set point is 60C, with a permitted tolerance of +/-5C (i.e. 55C to 65C), held for a duration of 4 hours (240 minutes).
3. Equipment qualification. Each oven must undergo preventive maintenance and recalibration on a 12-month schedule. A batch may not be released if it was processed on equipment whose most recent service is more than 12 months prior to the process date; such batches must be flagged for deviation review.
4. Temperature excursion handling. Any excursion above 65C for any duration must be logged as a deviation. Moisture content testing (spec: 0.0-0.5%) is the primary release test sensitive to over-drying caused by a high-temperature excursion.
5. This is a synthetic document generated for the DevGuard prototype and does not describe a real product or process.